



Naturally Radiant,
Eco-Friendly Rejuvenation!



Email: info@lpgclinicswholesale.com
Web: www.lpgclinicswholesale.com

iLLUMA PLLA

200mg x 1vial | 500mg x 1vial



poly l-lactid

iLLUMA PLLA rejuvenates the skin by stimulating collagen production, resulting in natural volume & firm, healthy skin

Composed of porous micro-particles

200mg x vial

500mg x vial



20~63µm

63~125µm

- Promotes tissue regeneration and excellent collagen production with minimal application
- Maintains a natural volumizing effect for a long time
- Consistent particle shape enables smooth injection without needle clogging

It stimulates collagen production within the face and body, thereby preserving the innate suppleness and elasticity of the skin



iLLUMA PLLA Product Features	
Volume :	200mg / 500mg
Dilution Volume :	8 ml / 20ml
Cannula :	iLLUMA PLLA 200 (Use with Cannula 27~28G) iLLUMA PLLA 500 (Use with Cannula 23~25G)
Main ingredients :	PLLA (Poly L-lactid) CMC Compound Amount
Number of treatment :	1-3 times
Duration :	18-24 month

Advantages of iLLUMA PLLA

The focus on natural and safe results has led to increased consumer interest, attracting and retaining patients.



Safety

The product, utilizing the NST method, is eco-friendly and solvent-free, ensuring safe and stable facial procedures.

Natural Outcomes

Made from natural corn-derived polymer, the product provides a voluminous effect and stimulates collagen regeneration.

Time-Efficiency

The significantly reduced dilution time cuts preparation to just 20 minutes, enhancing satisfaction and efficiency for both users and patients